BAMBANKER™ hRM

Serum-free, Xeno-free Cell Freezing Medium

<For Research Use Only>

BAMBANKER™ hRM is a xeno-free, serum-free cryopreservation medium used for cryopreserving variety of human derived cells including T lymphocytes. The cells are cryopreserved in a ultra low temperature freezer immediately after suspending with BAMBANKER™ hRM.

Features

- <1> BAMBANKER™ hRM is a ready-to-use medium, directly from a container without any preparation or dilution.
- <2> A controlled-rate freezer is not required.
- <3> BAMBANKER™ hRM allows rapid freezing and preserving cells immediately after placing the cells suspended with BAMBANKER™ hRM in a ultra low temperature freezer at -80°C.
- <4> BAMBANKER™ hRM is xeno-free, contains no serum.

An example of standard cryopreservation procedure

- 1) Pellet the cells in logarithmic growth phase.
 - 2) Suspend the collected cells with BAMBANKER™ hRM at a concentration of $5 \times 10^5 \sim 1 \times 10^7$ viable cells/mL and dispense aliquots of this suspension into cryogenic tubes
Then, freeze and preserve the cells in cryogenic tubes at -80°C without preliminary freezing.
Transfer the frozen cells into liquid nitrogen as required.
 - 3) In use, thaw the frozen cells in a cryogenic tube quickly in a constant-temperature chamber or water bath.
 - 4) Suspend the thawed cells into approximately 10 times the volume of fresh culture medium.
Mix well and spun down at 1,200rpm for 3 minutes, remove the supernatant.
Now the cells are ready for culturing and following procedures.
- ★ It is quite important to collect cells growing at logarithmic phase in order to secure a high rate of viable cells after thawing.

Sterility test

- Endotoxin: Colorimetric method
- Mycoplasma: Culture method
- Fungi and bacteria: membrane filtration method

Applications

Human T cells, Human iPS cells, Common Marmoset ES cells, Human iPS-derived livercells, Human iPS-derived neural cells, Human iPS-derived cardiomyocytes etc.

Precautions

- ◆ Use only for research purposes, not for human uses.
- ◆ Prior to use the product, perform a confirmation test of freezing and preservation on the cells under study.
- ◆ The manufacturer shall not be responsible for any accident of damage caused by the use of this product.
- ◆ For further assistance on the use of this medium, please contact the local distributor.

Genetics Catalog No.	Description	Package Size	Storage	Expiration period
BBH01	BAMBANKER™ hRM Serum-free Cell Freezing Medium, containing 10% DMSO	20mL	Keep at 2~10°C	Two years from manufacturing
BBH04		100mL		

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