

# BAMBANKER™ hRM

## Serum-free, Xeno-free Cell Freezing Medium

<For Research Use Only>

**BAMBANKER™ hRM** is a xeno-free, serum-free cryopreservation medium used for cryopreserving variety of human derived cells including T lymphocytes. The cells are cryopreserved in a ultra low temperature freezer immediately after suspending with **BAMBANKER™ hRM**.

### Features

- <1> **BAMBANKER™ hRM** is a ready-to-use medium, directly from a container without any preparation or dilution.
- <2> A controlled-rate freezer is not required.
- <3> **BAMBANKER™ hRM** allows rapid freezing and preserving cells immediately after placing the cells suspended with **BAMBANKER™ hRM** in a ultra low temperature freezer at -80°C.
- <4> **BAMBANKER™ hRM** is xeno-free, contains no serum.

### An example of standard cryopreservation procedure

- 1) Pellet the cells in logarithmic growth phase.
  - 2) Suspend the collected cells with **BAMBANKER™ hRM** at a concentration of  $5 \times 10^5 \sim 1 \times 10^7$  viable cells/mL and dispense aliquots of this suspension into cryogenic tubes  
Then, freeze and preserve the cells in cryogenic tubes at -80°C without preliminary freezing.  
Transfer the frozen cells into liquid nitrogen as required.
  - 3) In use, thaw the frozen cells in a cryogenic tube quickly in a constant-temperature chamber or water bath.
  - 4) Suspend the thawed cells into approximately 10 times the volume of fresh culture medium.  
Mix well and spun down at 1,200rpm for 3 minutes, remove the supernatant.  
Now the cells are ready for culturing and following procedures.
- ★ It is quite important to collect cells growing at logarithmic phase in order to secure a high rate of viable cells after thawing.

### Sterility test

Endotoxin: Colorimetric method  
Mycoplasma: Culture method  
Fungi and bacteria: Agar streak test

### Applications

Human T cells, Human iPS cells, Common Marmoset ES cells, Human iPS-derived livercells,  
Human iPS-derived neural cells, Human iPS-derived cardiomyocytes etc.

### Precautions

- ◆ Use only for research purposes, not for human uses.
- ◆ Prior to use the product, perform a confirmation test of freezing and preservation on the cells under study.
- ◆ The manufacturer shall not be responsible for any accident of damage caused by the use of this product.
- ◆ For further assistance on the use of this medium, please contact the local distributor.

Genetics Catalog No.	Description	Package Size	Storage	Expiration period
BBH01	<b>BAMBANKER™ hRM</b> Serum-free Cell Freezing Medium, containing 10% DMSO	20mL	Keep at 2~10°C	Two years from manufacturing
BBH03		120mL		

### Manufactured by:

**GC LYMPHOTEC Inc.**  
18-4, Fuyuki, Koto-Ku, Tokyo, 135-0041, Japan  
Tel: 81-3-3630-2530/ Fax: 81-3-3630-2531  
E-mail: info@lymphotec.co.jp  
http://www.lymphotec.co.jp

### Distributed by:

**NIPPON Genetics Co., Ltd**  
Mariaweiler Strasse 28-30, 52349 Dueren Germany  
Tel: 49-2421-55496-0/ Fax: 49-2421-55496-11  
E-mail: info@nippongenetics.eu  
http://www.nippongenetics.eu