# BAMBANKER™

## Serum-free type cell freezing medium

<For Research Use Only>

**BAMBANKER™** is a serum-free medium for long-term freezing at -80°C and preservation of valuable culture cells such as tumor cells and normal cells.

### Features

- <1> Ready-to-use medium for preservation of cells, <2> Usable without diluting, <3> No program freezer is required, <4> Rapid and long-term freezing and preservation in a deep freezer (-80°C), <5> Serum-free

### Operating Procedure

- 1) Centrifuge the cells in the logarithmic growth phase, and collect  $5 \times 10^5 \sim 1 \times 10^7$  cells.
- 2) Suspend the collected cells in 1 mL of this medium, and place the cells in cryotubes for freezing and preservation. Then, freeze and preserve the cells at -80°C without preliminary freezing. The frozen cells at -80°C can be preserved subsequently in liquid nitrogen.
- 3) Thawing must be carried out quickly in a constant-temperature chamber or bath.
  - ★ The cells must be frozen in the logarithmic growth phase.

### Sterility test

Endotoxin: Colorimetric method  
Mycoplasma: Culture method  
Fungi and bacteria: Agar streak test

### Applications <Cells preserved well by cryopreservation test>:

- P3U1 (mouse myeloma cell line), K562 (human leukemia cell line), human gastric epithelial cells, human  $\gamma\delta$ T cells
- Daudi (human B cell line), PC12 (rat-derived adrenal pheochromocytoma), human B cell line
- OKT4 (mouse hybridoma), monkey B cell line, Activated lymphocyte derived from human peripheral blood, Activated lymphocyte derived from mouse spleen

### Precautions

- ◆ Use only for research purposes, not for human use.
- ◆ Before using this medium, perform the confirmation test on the cells under study.
- ◆ The manufacturer shall not be responsible for any accident or damage caused by the use of this product.
- ◆ If you have any questions about the use of this medium, please contact the local distributor.

| Description                                                                  | Package Size | Storage           | Expiration period               |
|------------------------------------------------------------------------------|--------------|-------------------|---------------------------------|
| <b>BAMBANKER™</b><br>Serum-free Cell Freezing<br>Medium, containing 10% DMSO | 120mL        | Keep at<br>2~10°C | Two years from<br>manufacturing |
|                                                                              | 20mL × 5     |                   |                                 |

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